

# Ion milling of a 100 nm-thick Silicon Nitride membrane

Briosne Frejaville Clemence  
Mau Adrien

Info for all following slides:

- FIB miling time  $< 1$  s
- SEM images were made without any coating

**Main result:** Resizing of holes down to  $\sim 10$  nm diameter is possible (challenge is to measure)

# Smallest holes (10 nm is resolution spec for IonLine)




Figure 1 :

Smallest Ga-ion milled holes.  
Holes with diameter around 20nm.




Figure 2 :

Smallest Ga-ion milled holes.  
Holes with diameter around 20nm.

# Back-side view




Figure 3: Back-side view of the hole. The 40nm wide hole (front view) is 20 nm (back view).

# Resizing of holes: easy!




Figure 4: One of the two holes of 80 nm was closed down to 20 nm hole. The shape of the final hole is also more circular.